









CondenSyn MAXI for reflux or for distillation

Details of the different configurations available



CondenSyn MAXI Reflux Inner

CondenSyn MAXI Outer

CondenSyn MAXI Distillation Inner



The image (left) shows the CondenSyn MAXI with the distillation set up, as it would be used in conjunction with a lab reactor vessel.

The outer component is identical to that in the reflux setup, but the inner component has a small drain arm at the bottom that has a GL 18 thread on it.



- GB-C-MAXI-350-P-B34 is a 350mm length in reflux configuration.
- GB-C-MAXI-350-S-B34 is a 350mm length in distillation configuration.
- GB-C-MAXI-450-P-B34 is a 450mm length in reflux configuration.
- •MAXI-450-S-B34 is a 450mm length in distillation configuration.

Contact us now for further details & pricing



42 à 48 Bd de Polangis - BP 260 - 94502 Champigny-sur-Marne Cedex Tél.: 01 48 83 21 76 - Fax.: 01 48 83 51 01 info@cloup.fr www.cloup.fr